

# How to optimize your chromatography

## Factors influencing Speed and Resolution

**Author:**  
**Michael Frank**  
Application Chemist  
Agilent Technologies  
Waldbronn, Germany


**Speaker:**  
**Rob Solazzo**  
Asia Pacific Industry Marketing  
Manager  
Pharma/Biotech

## Factors to improve

- Speed
  - Run-time
  - Cycle-time
- Resolution
- Costs
  - Labor costs
  - Instrument costs
  - Solvent consumption
  - Waste production

## The Effect of Particle Size

Van-Deemter curve for different particle sizes









## The Effect of Particle Size

Reducing the particle size will yield a...

- à decreased  $H_{min}$  (better efficiency)
- à at higher  $u$  (higher speed)
- à with decreased slope (less loss of efficiency when increasing  $u$ )

## The Effect of Particle Size

*Development of stationary phases*

Year	Particle size	Plates/50 mm
1969	 100 $\mu\text{m}$	170
1973	 57 $\mu\text{m}$ (pellicular)	350
1975	 10 $\mu\text{m}$	2,000
1985	 5.0 $\mu\text{m}$	4,000
1992	 3.5 $\mu\text{m}$	7,500
2003	 < 1.8 $\mu\text{m}$	12,000

## The Effect of Particle Size

*The big drawback of using smaller particles:*


Pressure increases with the square of  $1/d_p$

$$\Delta p \propto \frac{4 \cdot F \cdot L_c \cdot h}{pd_c^2 \cdot d_p^2}$$

$F$  = flow,  $L_c$  = column length,  $h$  = viscosity,  $d_c$  = column diameter

## The Effect of Particle Size

Pressure dependence of length, diameter and particle size:



Flow = 1 mL/min, viscosity = 0.75 cP


Page 7

Agilent Technologies

Group/Presentation Title  
Agilent Restricted  
Month ##, 200X

## The Effect of Particle Size

Tricks to overcome this – optimize the particle size distribution (PSD)




Page 8

Agilent Technologies

Group/Presentation Title  
Agilent Restricted  
Month ##, 200X

## The Effect of Particle Size

Tricks to overcome this – optimize the particle size distribution (PSD)



Page 9

Agilent Technologies

Group/Presentation Title  
 Agilent Restricted  
 Month ##, 200X

## How to use smaller particle sizes to improve existing method?

- *Use a shorter column:*  
 Compounds will **elute earlier** but with **preserved resolution** (higher efficiency  $N$  of smaller particles).  
 Also, the linear velocity can be increased more than with larger particles because lower loss of efficiency with increasing velocity of smaller particles.

- *Keep column length constant:*  
**Resolution  $R_S$  will increase** with same retention time.

$$N \propto \frac{1}{d_p} \quad R_S \propto \sqrt{N}$$

Page 10

Agilent Technologies

Group/Presentation Title  
 Agilent Restricted  
 Month ##, 200X

***How to increase the speed of a conventional method?***

## Increasing speed

- *Length scaling – to preserve efficiency adjust the particle size*

$$N = \frac{L_{c1}}{d_{p1}} = \frac{L_{c2}}{d_{p2}} \quad \Rightarrow \quad L_{c2} = \frac{L_{c1} \cdot d_{p2}}{d_{p1}}$$

Example – present method running on a 250 x 4.6mm, 5.0µm:

$$L_{c2} = \frac{250mm \cdot 1.8\mu}{5\mu} = 90mm$$

So, use the closest available size – a 100 x 4.6mm, 1.8µm column

## Increasing speed

- *Length scaling – adjust the time*

$$\frac{L_{c2}}{L_{c1}} = \frac{t_2}{t_1} \quad \Rightarrow \quad t_2 = \frac{L_{c2} \cdot t_1}{L_{c1}}$$

Example – scaling a 30min-method from a 250 mm to 100 mm column:

$$t_2 = \frac{100\text{mm} \cdot 30\text{min}}{250\text{mm}} = 12\text{min}$$

## Increasing speed

- *If wanted, scale the diameter, too – but, adjust the flow:*

$$\left(\frac{d_{c2}}{d_{c1}}\right)^2 = \frac{F_2}{F_1} \quad \Rightarrow \quad F_2 = \left(\frac{d_{c2}}{d_{c1}}\right)^2 \cdot F_1$$

**Note:** With a given max. flow of a LC system a higher linear velocity can only be achieved with a smaller column ID!

Example – scaling a 1.4ml/min method from 4.6 mmID to 3 mmID:

$$F_2 = \left(\frac{3.0\text{mm}}{4.6\text{mm}}\right)^2 \cdot 1.4\text{ml/min} = 0.595\text{ml/min}$$

## Increasing speed

- *Optimizing the linear velocity (flow)*

Smaller particles have an increased  $u_{opt}$ ; so, flow should be increased

$$\frac{F_1}{F_2} = \frac{t_2}{t_1} = \text{const.} \quad \Rightarrow \quad t_2 = \frac{F_1 \cdot t_1}{F_2}$$

Example – increasing flow from 1.4ml/min to 3.6ml/min:

$$t_2 = \frac{1.4 \text{ ml/min} \cdot 12 \text{ min}}{3.6 \text{ ml/min}} = 4.7 \text{ min}$$

Note: if the linear velocity is increased beyond the optimal velocity a reduction in efficiency must be expected!

## Increasing speed

- *Optimizing the linear velocity (flow) – cont.*

For gradient separations the gradient slope  $s$  should be preserved

$$s = \frac{\Delta\%B}{(F_1 \cdot t_{g1}) / V_{M1}} = \frac{\Delta\%B}{(F_2 \cdot t_{g2}) / V_{M2}}$$

This describes the gradient slope as  $D\%B$  per number of column void volumes ( $V_M$ ).



## Increasing speed

- Finally, adjusting the injection volume

To achieve same peak heights (LOD, LOQ) the injection volume needs to be scaled according to the column void volumes

$$\frac{V_{inj2}}{V_{inj1}} = \frac{V_{M2}}{V_{M1}}$$


Assuming an identical injection volume, with larger columns, i.e. with a bigger void volume, the peaks would dilute more, resulting in a smaller peak height and a lower detection limit.



Same sample amount but different height if the injection volume is not adjusted!

## Increasing speed

- Summary



And of course:  $N = \frac{L_{c1}}{d_{p1}} = \frac{L_{c2}}{d_{p2}} = \text{const.} \quad !!!$

$$F_2 \leq F_{\max} = \frac{p_{\max} d_{c2}^2 \cdot d_{p2}^2 \cdot p_{\max}}{4 \cdot \Phi \cdot L_{c2} \cdot h}$$


Max. flow  $F_{\max}$  dependent on max. system pressure  $p_{\max}$ , Column ID  $d_c$  & Length  $L_c$ , Particle size  $d_p$ , Viscosity  $h$  and the flow resistance factor  $F$  of the column

**How to increase the resolution of a conventional method?**

## Increasing Resolution

- The resolution equation:

$$R_s = \frac{\sqrt{N}}{4} \cdot \left[ \frac{a-1}{a} \right] \cdot \left[ \frac{k_1}{k_2+1} \right]$$




## Increasing Resolution

**Selectivity (a)** helps best but:

- Is difficult to predict (experience helps, model retention)
- Requires, software supported, optimization for separation of multi-component mixtures (*ChromSword*, *DryLab*)

### Example:

Same sample, analyzed with different conditions but always same temperature and 5-95%B gradient in the same time.



Page 21

Agilent Technologies

Group/Presentation Title  
Agilent Restricted  
Month ##, 200X

## Increasing Resolution


High **plate number (N)** by decreasing the particle size of a column provides:

- Sharp and narrow peaks
- Better detection (better S/N)
- Peak capacity to resolve complex samples

But, resolution increases only with the **square root of the plate number** and plate number increase is **limited by experimental conditions** (analysis time, pressure)!

### Example:

Same sample, analyzed on a 4.6mm x 50mm column but with different  $d_p$ .



Page 22

Agilent Technologies

Group/Presentation Title  
Agilent Restricted  
Month ##, 200X

## Increasing Resolution

Increasing the **retention factor ( $k$ )**

is very helpful between 5 and 10 but beyond that only off insignificant improvement.

In gradient separations  $k$  can be influenced by the gradient slope, i.e. flow and/or gradient time.

$$k^* = \frac{t_g \cdot F}{1.15 \cdot S \cdot \Delta\% B \cdot V_M} \quad (\text{average } k \text{ for gradient separation})$$

So, by just increasing the flow in a gradient separation  $k$  becomes bigger and the resolution will improve!

## All scaling calculations to transfer methods to sub-2 $\mu$ m particles are done by the Agilent MethodTranslator

The screenshot displays several overlapping windows from the Agilent MethodTranslator software. The main window in the foreground is titled 'Agilent MethodTranslator' and shows a 'Method Translator' dialog box. It includes fields for 'Original Method' (Agilent 1100 Series LC) and 'Target Method' (Agilent 1200 Series HPLC). The 'Scaling Factor' is set to 1.0. Below these fields are two data tables. The first table, 'Original Method', lists parameters like 'Flow Rate (mL/min)', 'Injection Volume (μL)', and 'Retention Time (min)'. The second table, 'Target Method', lists 'Flow Rate (mL/min)', 'Injection Volume (μL)', and 'Retention Time (min)'. To the right, there is a 'Method Translator' window with a 'Scaling Factor' of 1.0 and a 'Method Translator' button. In the bottom right corner, there is a red text box that reads 'For free - P/N 5989-5130EN'.

## Features of the Agilent MethodTranslator

Basic mode with certain pre-set parameters:

Enter the parameters of your existing method and the parameters of the desired column you would like to convert to.

Page 25

Agilent Technologies

Group/Presentation Title  
Agilent Restricted  
Month ##, 200X

## Features of the Agilent MethodTranslator

Advanced mode – all calculation parameters in your hands:

More to enter but much more information returned.

Page 26

Agilent Technologies

Group/Presentation Title  
Agilent Restricted  
Month ##, 200X

# Features of the Agilent MethodTranslator

## Viscosity Tables

Temperature Dependent Viscosity [cP] of Water / Acetonitrile Mixtures


Temp [°C]	0	5	10	15	20	25	30	35	40	45	50	55	60	65	70	75	80	85	90	95	100	
20	1.00	1.02	1.12	1.18	1.08	1.05	1.03	1.01	0.99	0.97	0.95	0.93	0.91	0.89	0.87	0.85	0.83	0.81	0.79	0.77	0.75	0.73
25	0.98	1.01	1.10	1.16	1.06	1.03	1.01	0.99	0.97	0.95	0.93	0.91	0.89	0.87	0.85	0.83	0.81	0.79	0.77	0.75	0.73	0.71
30	0.97	1.00	1.09	1.15	1.05	1.02	1.00	0.98	0.96	0.94	0.92	0.90	0.88	0.86	0.84	0.82	0.80	0.78	0.76	0.74	0.72	0.70
35	0.96	0.99	1.08	1.14	1.04	1.01	0.99	0.97	0.95	0.93	0.91	0.89	0.87	0.85	0.83	0.81	0.79	0.77	0.75	0.73	0.71	0.69
40	0.95	0.98	1.07	1.13	1.03	1.00	0.98	0.96	0.94	0.92	0.90	0.88	0.86	0.84	0.82	0.80	0.78	0.76	0.74	0.72	0.70	0.68
45	0.94	0.97	1.06	1.12	1.02	0.99	0.97	0.95	0.93	0.91	0.89	0.87	0.85	0.83	0.81	0.79	0.77	0.75	0.73	0.71	0.69	0.67
50	0.93	0.96	1.05	1.11	1.01	0.98	0.96	0.94	0.92	0.90	0.88	0.86	0.84	0.82	0.80	0.78	0.76	0.74	0.72	0.70	0.68	0.66
55	0.92	0.95	1.04	1.10	1.00	0.97	0.95	0.93	0.91	0.89	0.87	0.85	0.83	0.81	0.79	0.77	0.75	0.73	0.71	0.69	0.67	0.65
60	0.91	0.94	1.03	1.09	0.99	0.96	0.94	0.92	0.90	0.88	0.86	0.84	0.82	0.80	0.78	0.76	0.74	0.72	0.70	0.68	0.66	0.64

Temperature and composition dependent for H<sub>2</sub>O/ACN and H<sub>2</sub>O/MeOH

*Does it work?*

## Does it work? - Example

Analysis of impurities of an active pharmaceutical ingredient by conventional HPLC (4.6mmID x 250mm, 5.0µm):



Page 29

Agilent Technologies

Group/Presentation Title  
Agilent Restricted  
Month ##, 200X

## Does it work?

Converting to a 4.6 mmID x 100 mm, RRHT column:

Method Information for Agilent 1100 Series LC:

- System: Agilent 1100 Series LC
- Column: 4.6 mmID x 250 mm
- Flow Rate: 1.40 mL/min
- Injection Vol: 10.00 µL
- Pressure: 24.77

Method Information for Agilent 1200 Series RRLC:


- System: Agilent 1200 Series RRLC
- Column: 4.6 mmID x 100 mm
- Flow Rate: 2.500 mL/min
- Injection Vol: 10.00 µL
- Pressure: 24.000

Page 30

Agilent Technologies

Group/Presentation Title  
Agilent Restricted  
Month ##, 200X

## Does it work?



4.6 mmID x 100 mm, 1.8µm Zorbax SB C18

Speed Optimized


0.00 min	5% B
4.33 min	90% B
4.98 min	90% B
4.99 min	5% B
6.5 min	5% B

Page 31

Agilent Technologies

Group/Presentation Title  
Agilent Restricted  
Month ##, 200X

## Does it work?



Conventional HPLC

Resolution optimized 1.8µm


Page 32

Agilent Technologies


Group/Presentation Title  
Agilent Restricted  
Month ##, 200X




## Does it work? - Resolution




## Does it work? – Run Time



## Does it work? – Solvent Consumption



## Does it work? - Pressure Prediction of the MethodTranslator\*



\*) calculated with the standard setting of the basic mode. Note: at very high flow rates larger deviations can be expected as the backpressure of capillaries are not taken into account.

## The MethodTranslator - What about reality?

- **You have to take additional backpressure of capillaries into account.**  
*Smaller particles  $\Rightarrow$  higher Efficiencies  $\Rightarrow$  higher demands on dead volume, BUT: backpressure induced by the column is proportional to  $d_p^2$  but induced by capillaries proportional to  $d_c^4$*
- **Real columns have a distribution of the particle size.**  
*We have a higher amount of larger particles around the mean of  $1.8\mu\text{m}$ , other manufactures have a higher amount of fines around their mean particle size.*
- **What is the real temperature of the column?**  
*You have to consider energy lost by transmission but energy gain from friction! So, what is the true temperature of the column, the true viscosity?*
- **Don't forget about the real measured Efficiency!**  
*The Efficiency of the column is deminished (sometimes dramatically) by extracolumn volumes which cause bandspreading, injection volumes (not always an ideal solute can be used) and by the detection (flow cell, acquisition rate). Depending on original method and new method these contributions can be different!*



## What is the MethodTranslator and what is it not?

***The MethodTranslator is a tool to help you in converting established methods to smaller particle size columns and utilize their advantages. It gives you a fairly good starting point, but...***

***...the MethodTranslator is **not** a method development software like ChromSword or DryLab!***



## Summary

- **Sub-2 $\mu$ m particle columns allow to get either shorter analysis times and/or higher resolutions**
- **Time and solvent savings can be significant**
- **Transfer is very easy – the Agilent MethodTranslator is a tool to support you**
- **Nowadays, methods using sub-2 $\mu$ m columns are as robust as conventional HPLC methods**

For more information visit:  
[www.agilent.com/chem/1200rr](http://www.agilent.com/chem/1200rr)

Agilent Technologies

### 1200 系列快速

Agilent 1200 系列快速 HPLC 系统，为您提供更短的分析时间，更高的分辨率，更少的溶剂消耗，以及更低的运行成本。该系统采用亚 2 微米颗粒填料，具有卓越的性能，是您实验室的理想选择。

产品特点

- 亚 2 微米颗粒填料，提供更短的分析时间，更高的分辨率，更少的溶剂消耗，以及更低的运行成本。
- 优化的流路设计，减少死体积，提高峰容量。
- 兼容多种溶剂，满足不同应用需求。
- 易于操作和维护，降低运行成本。

联系我们

Agilent (CN) 客服热线  
400-887-7237


## Backup

Method	ID [mm]	Length [mm]	Particle Size [µm]	Flow [mL/min]	Run Time [min]	Solvent consumption [mL]	Solvent cons. vs. Conv.	Rs (A, D)	Rs (Main, A)	Rs (C, Main)	Rs (B, C)	predicted pressure *) [bar]	measured pressure [bar]	Deviation of predicted pressure
Conventional HPLC	4.6	250	5.0	1.400	23.00	32.20	100%	19.32	2.49	18.50	4.19	94.77	103.50	8.43%
4.6 Simple conversion	4.6	100	1.8	1.400	9.20	12.88	40%	19.83	2.30	16.88	4.32	292.50	309.00	5.34%
4.6 Speed optimized	4.6	100	1.8	2.585	4.98	12.87	40%	17.29	2.12	15.99	3.88	540.00	548.00	1.46%
4.6 Speed optimized 80 ° C **)	4.6	100	1.8	4.500	2.86	12.87	40%	14.53	1.19	12.50	3.41	540.00	569.00	5.10%
4.6 Resolution optimized	4.6	100	1.8	2.585	9.20	23.78	74%	23.37	2.83	20.22	5.08	540.00	549.00	1.64%
5.0 Speed optimized	3.0	100	1.8	1.100	4.98	5.48	17%	16.77	2.12	16.13	3.76	540.00	547.00	1.28%
2.1 Speed optimized	2.1	100	1.8	0.539	4.98	2.68	8%	14.62	2.07	14.83	3.28	540.00	514.00	5.06%
2.1 Speed optimized **)	2.1	50	1.8	1.077	1.25	1.35	4%	7.18	1.09	7.71	1.43	540.00	515.00	4.85%

\*) calculated with the standard setting of the basic mode. Note: at very high flow rates larger deviations can be expected as the backpressure of capillaries are not taken into account.


\*\*\*) Note: the resolution for the critical pair is below the usually accepted limit of 2 for such an analysis.

## Backup



Comparison of experimental kinetic plots for high temperature LC (5µm, 30 and 150 °C) and high pressure UPLC (1.7µm, 40 °C). Maximum pressures: 400 and 1000 bar (UPLC), 400 bar (LC). Values for  $\phi^0$ , A-C,  $\eta$  and  $D^0$  are listed in Table 5. Columns and mobile phases: Blaze C18 200 (150mm×4.6 mm, 5µm with 40/60 acetonitrile/water and Acquity BEH C18 (100mm×2.1 mm, 1.7µm) with 30/70 acetonitrile/water.

## Backup



Theoretical plate-height curves for 5m particles, illustrating the effect of temperature on plate height and linear velocity. The plots were calculated with Eq. (12) assuming values of:  $A = 0.66$ ,  $B = 3.00$  and  $C = 0.05$ . Diffusion coefficients for phenol in 40/60 acetonitrile/water were calculated at each temperature according to Wilke-Chang and Eqs. (11) and (13) [30,19]:  $1 \times 10^{-9}$  m<sup>2</sup>/s (25 °C),  $2.4 \times 10^{-9}$  m<sup>2</sup>/s (75 °C),  $4.5 \times 10^{-9}$  m<sup>2</sup>/s (125 °C),  $7.7 \times 10^{-9}$  m<sup>2</sup>/s (175 °C).